

INSTALLING CAE MAESTRO STANDALONE FOR MAC


This guide describes how to install CAE Maestro standalone software on any personal Mac computer. The software can be downloaded from any Mac computer with an Internet connection.

Note: To install updates for the Instructor Workstation, see the corresponding *CAE Maestro Software Conversion* guide.

Installation Steps	
1	Download and Install CAE Maestro Standalone
2	Activate the Software License and Launch CAE Maestro Standalone


Step 1: Download and Install CAE Maestro Standalone

- Locate **cae-maestro-standalone-darwin-x.x.x.x.pkg** in the **Downloads** folder and double-click on it.




Downloads Folder

- b. The Maestro Installer window will appear. Click **Continue**.




Maestro Installer Window

- c. When the Important Information window appears, click **Continue**.




Maestro Installer Window

d. The Software License Agreement window will appear. Click **Continue**.




Software License Agreement Window - Click Continue

e. To continue installing the software, click **Agree**.




Software License Agreement Window - Click Agree

- f. The Destination Drive window will appear. Click **Continue**.




Destination Drive Window

- g. The installation window will appear. Click **Install**.




Installation Window

- h. Enter the **username** and **password** of a computer administrator, then click **Install Software**.




Administrator Authorization to Install Software

Wait while the CAE Maestro Standalone software installs.




CAE Maestro Software Installing

- i. A CAE Maestro Standalone shortcut will appear on the Desktop. Click **Close** to close the installation window.



CAE Maestro Software Installation Successful


- j. A prompt may appear. Click **Move to Trash**. Otherwise, drag the Installer to the Trash.



Prompt - Move Installer to Trash

The Installer Window closes and the cae-maestro-standalone-darwin-x.x.x.x.pkg file is moved to the trash. Close the Downloads folder by clicking on the close button (red circle) in the upper-left corner of the Downloads folder window.

Close
button



Close Button


k. Empty the Trash by selecting **Finder > Empty Trash**.

Finder
menu



Finder Menu - Empty Trash

- l. A prompt will appear. Click **Empty Trash**.



Prompt - Empty Trash


- m. Click the **Apple** icon in the upper-left corner of the Desktop, then select **Restart**.

Apple icon



Apple Menu - Restart

n. A prompt will appear. Click **Restart**.



Prompt - Restart

Step 2: Activate the Software License


Each patient simulator includes one free license for CAE Maestro Standalone. You must activate the software license either online or by calling CAE Customer Service before the update takes effect.

Prior to activating the software license in Maestro, you must locate the license key included in the simulator documentation or provided by CAE via email.

Note: If you cannot access our server via Internet, the computer will generate a 25-character Simulator Number to be used when contacting Customer Service.


To launch the Maestro Standalone software:

- a. From the desktop, click on the new **CAE Maestro Standalone** shortcut icon to open the new version of CAE Maestro Standalone.




CAE Maestro Standalone Shortcut Icon

- b. From the CAE Maestro Home screen, click the **CAE Maestro** icon.




Maestro Home Screen

- c. Select the desired simulator in the dropdown box and click **Continue**.



Select Your Simulator Screen

The Simulator Software Unlicensed message appears.




Simulator Software Unlicensed Message

- d. Click **License Manager** to activate the license online or offline.

Online Activation


To activate the software license online:

- a. From the License Manager screen, on the Simulation Software tab, click **Activate**.



License Manager Screen

The Activate Product screen appears.




Activate Product Screen

- b. Enter the license key provided in your simulator documentation or email from CAE and click **Activate**.

The Success message appears. Close the Success window.

- c. Confirm the software license is activated on the License Manager screen. A “Valid License” indicator appears in line with the software.




Activated License

Note: Follow the same procedure to license the Maestro Physiology.

Offline Activation

To activate offline:


- a. From the License Manager screen, on the Simulation Software tab, click **Activate**.



License Manager Screen


The Activate Product Offline screen appears.

- b. Call Customer Service (see *Customer Service and Training on page 10*) with your License Key and 25-character Simulator Number.
You will receive an activation key.




Activate Product Offline Screen

- c. Enter the license key provided. The Incorrect format message will appear until the full license key is entered.



Activate Product Offline Screen


- d. When completed, click **Activate**.



Activate Product Offline Screen

The Success message appears.

- e. Confirm the software license is activated on the License Manager screen. A “Valid License” indicator appears in line with the software.



Activated License

Note: Follow the same procedure to license the CAE Maestro Physiology.

- f. To launch CAE Maestro Standalone, on the Desktop double-click the **CAE Maestro Standalone** shortcut.



CAE Maestro Standalone Shortcut

Customer Service and Training

For assistance, contact customer service.

Customer Service Headquarters - United States

Monday - Friday from 7:00 a.m. to 6:00 p.m. ET

Toll Free: +1 (866) 462-7920

Fax: +1 (941) 342-5600

Email: srqcustomerservice@cae.com

Website: www.caehealthcare.com

Customer Service - Canada

Monday - Friday from 8:00 a.m. to 5:00 p.m. ET

Toll Free: +1 (877) 223-6273

Email: can.service@cae.com

Customer Service - Europe, Middle East and Africa

Monday - Friday from 8:00 a.m. to 5:00 p.m. CET

Phone: +49-(0) 6131 4950354

Fax: +49 (0) 6131 4950351

Email: international.service@cae.com

Customer Service - United Kingdom and Ireland

Monday - Friday from 9:00 a.m. to 5:00 p.m. GMT

Phone: +44(0)800-917-1851

Email: uk.service@cae.com

Customer Service - Latin America

Monday - Friday from 9:00 a.m. to 5:00 p.m. BRT / BRST

Phone: +55 11 5069-1510

Email: la.service@cae.com

Customer Service - Asia Pacific

Monday - Friday from 8:00 a.m. to 5:00 p.m. CET

Phone: +49 (0) 6131 4950354

Email: ap.service@cae.com

